

Bed Exit Mat

Part No. **QT633**



Overview

- Designed to trigger a call by sensing the pressure change when a patient exits a bed.
- Connects to Quantec via the remote jack socket on a QT637 Enuresis/General Purpose Interface Socket or QT636 Interface Unit.
- To reduce the risk of nuisance calls, the QT637 includes circuitry that prevents a momentary change in pressure such as a patient quickly changing position from triggering a call.
- Typically positioned across and over a mattress, directly underneath the bed sheet.
- Measures 80cm x 10cm and is supplied with a 2m lead and normally closed right angled jack plug.

More Information

Technical Specifications

Compatibility	Must be used in conjunction with a QT637 Quantec Enuresis Sensor/General Purpose Interface Socket or QT636 Interface Socket.
Connections	6.35mm (quarter inch) jack plug.
Product dimensions (mm)	766mm W x 102mm H. Includes a 2m lead.
Construction & finish	Plastic.
IP Rating	TBC.
Weight	240g
Operating conditions/temperature	-5°C to +40°C. Max relative humidity: 95%.